Baxter

6208750ZA01

EN

TRAUMASTEM Oxidised Cellulose

Sterile resorbable haemostatic reticulum

Manufacturer:	BIOSTER, a.s., Tejny 621, 664 71 Veverská Bítýška, Czech Republic
Distributor:	Baxter Healthcare SA, 8010 Zurich, Switzerland
Composition:	Oxidized cellulose (100% active substance)
Indication group:	Haemostats (local)
Code UMDNS:	17 944

Characteristics:
Traumastem is a sterile, resorbable, biocompatible, haemostatic reticulum designed to stop capillary, venous and minor arterial bleeding. Traumastem provides fast and effective local haemostasis within 1 to 2 minutes depending on the type of surgery and the intensity of bleeding. High absorption ability as well as large surface area of oxidized natural cellulose contributes greatly to the process of haemostasis. After contact with blood oxidized cellulose changes into an adhesive gel which protects thrombocytes from being washed out of the wound, thereby stopping bleeding faster. It minimizes blood loss and complications perioperatively and in the early post-operative stages, eliminates the risk of blood transfer, has significant antibacterial activity, contributes to the acceleration of biochemical processes and facilitates the healing process. pH value of aqueous leach of product is 2.7 to 3.5. It is atraumatic and fully biodegradable. It is fully absorbed by the organism within 48 to 72 hours (depending on wound characteristics) and thereby product remaining in the tissue and/or cells and possible post-op irritation does not occur. Traumastem is well tolerated by the body and does not cause any undesirable immune reactions or sensitization even when used repeatedly. When in contact with body fluids/blood it does not increase its volume, does not expand into surrounding areas and does not cause any pressure on surrounding tissues or organs. Therefore it is not necessary to remove Traumastem from the area where applied. Traumastem has antimicrobial activity *in vitro* against these microorganisms in compliance with USP 31, article 51: *Staphylococcus epidermidis* *Escherichia coli* *Pseudomonas aeruginosa* *Candida albicans* *Aspergillus niger* *Staphylococcus aureus* *Proteus vulgaris* *Enterococcus faecalis* *Klebsiela aerogenes* *Bacteroides fragilis* *Clostridium perfringens*

Indications:
Traumastem is used to stop capillary, venous, small arterial bleeding as well as in the prevention of bleeding in the early post-operative stage. Traumastem represents a unique product for use in all areas where stoppage of diffuse bleeding is required for instance from resected areas of organ parenchyma, muscles or circumscribed body cavities (for example small pelvis) without complications. Traumastem is indicated for use in general and digestive surgery, plastic surgery, vascular and thoracic surgery, neurosurgery, orthopedics, gynaecology, urology, stomatology and other similar surgical fields. Its wide range of application enables the product to be used in surgeries carried out traditionally or laparoscopically, in robotic surgery, for endoscopic procedures as well as for one-day surgeries.

Contraindications:
Traumastem is not designed to be used to stop large arterial bleeding.

Directions for Use:
Traumastem is supplied in a sterile form. After being taken out of package the haemostatic reticulum is applied directly onto the dried bleeding area. Bleeding is usually stopped within 1 or 2 minutes depending on the type of surgery and the intensity of bleeding. It acts locally. Traumastem can either be applied to the required area by gentle pressure, attaching it with a sitch, sewing it onto a tissue or an organ or it can be used as a form of tamponade. It is perfectly absorbed by the body and does not need to be removed from the place of application, thereby preventing the recurrence of bleeding, client traumatization or interruption of the healing process. Traumastem can be cut to the size of the bleeding area. During neurosurgical procedures according to the surgeon's experience it is up to the surgeon's discretion whether or not to leave the reticulum in the surgical field.

Warning:
The product is made for a single usage only. Product sterility is ensured only by undamaged packaging. Although Traumastem has antibacterial effect it is not intended as a substitute for systemically administered therapeutic or prophylactic antimicrobial agents. The product cannot be resterilized. Haemostatic activity of Traumastem is more efficacious when it is applied dry, therefore it should not be moistened with water or physiological solution before application.

Description:
Traumastem is a sterile, resorbable, haemostatic reticulum of white or yellowish colour. It is designed for outer as well as inner usage.

Packaging:	10 pieces		
Dimensions:	1,5 cm x 5 cm 7,5 cm x 5 cm 12,5 cm x 5 cm 35 cm x 5 cm 20 cm x 10 cm		It is possible to customize dimensions on request.

Storage:
Store in original packaging, in a dry place at temperatures ranging from 5 to 25 °C.

	Do not use if package is damaged.
	Manufacturer.
	Distributor.
	Consult Instructions for Use.

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